







Materials available for use with CEREC

	Block / Material	Mpa	Use: Molar (Full Coverage)	Use: Molar (Partial Coverage)	Use: Premolar/Canine	Use: Incisors	Use: Temps	www.pattersondental.com link
	3M ESPE Lava Ultimate	210		✓	✓	✓		https://goo.gl/iX3H1Y
	3M ESPE MZ100	150		✓	✓	✓		https://goo.gl/5exhAK
	CEREC Guide 2							https://goo.gl/Q2rkd9
	Dentsply Celtra Duo	210-370		✓	✓	✓		https://goo.gl/gzD6JA
	GC America CERASMART	240	✓	✓	✓	✓		https://goo.gl/TnE9um

Materials available for use with CEREC

	Ivoclar e.max CAD	500	✓	✓	✓	✓	https://goo.gl/cV3gtF
	Ivoclar e.max CAD Abutment	500	✓	✓	✓	✓	https://goo.gl/rd1yLR
	Ivoclar e.max ZIRCAD	1200	✓	✓	✓	✓	https://goo.gl/oAvs2x
	Ivoclar Empress	160		✓	✓	✓	https://goo.gl/dLPks6
	Ivoclar Empress Multi	160		✓	✓	✓	https://goo.gl/6rVM3K
	Ivoclar Telio CAD	130	✓	✓	✓	✓	https://goo.gl/osh6sH





Materials available for use with CEREC

	Kuraray Katana	760	✓	✓	✓	✓	✓	https://goo.gl/z7vmMN
	Shofu Block HC	190		✓	✓	✓		https://goo.gl/45G8cr
	Sirona Bloccs C	120-160		✓	✓	✓		https://goo.gl/3NmneL
	Sirona Bloccs C In	120-160		✓	✓	✓		https://goo.gl/br5GsS
	Sirona Blocks C PC	120-160		✓	✓	✓		https://goo.gl/v6nVzB
	Sirona CEREC Zirconia	>900	✓	✓	✓	✓		https://goo.gl/2eNVJa

Materials available for use with CEREC

	Sirona InCoris TZI C	>900	✓	✓	✓	✓	https://goo.gl/8D6NCq
	Sirona InCoris ZI meso	>900	✓	✓	✓	✓	https://goo.gl/YJs8Xq
	Vita CAD Temp	>80	✓	✓	✓	✓	https://goo.gl/of5Jfk
	Vita CAD Temp Multi	>80	✓	✓	✓	✓	https://goo.gl/Yvy2NX
	Vita Enamic	150-160	✓	✓	✓	✓	https://goo.gl/Zddf4L
	Vita Enamic Multi	150-160	✓	✓	✓	✓	https://goo.gl/PuJFBf

Materials available for use with CEREC

	Vita Enamic IS	150-160	✓	✓	✓	✓	https://goo.gl/AUq6a2
	Vita Mark II	120-160		✓	✓	✓	https://goo.gl/sf2otx
	Vita RealLife	120-160		✓	✓	✓	https://goo.gl/owCjAa
	Vita Triluxe	120-160		✓	✓	✓	https://goo.gl/Ubf2Jm
	Vita Triluxe Forte	120-160		✓	✓	✓	https://goo.gl/LncAou